

# CITO Cushion Crease CITO Cushion Crease Plus


High-quality creasing profile for the processing of corrugated board

Especially developed for flatbed and rotary dies, for diemakers and corrugated board converters

- Excellent bonding capability
- Outstanding height and wear resistance
- Simple optical quality control through colour-coding
- Available as self-adhesive EasyFix version on user-friendly CMR reels
- The self-adhesive EasyFix technology with Fingerlift makes the application even more user-friendly and efficient




## CITO Cushion Crease – creasing profile



CITO Cushion Crease is used on the creasing rule to support compression of the corrugated board and to optimise folding behaviour, particularly along the corrugations.


### Cushion Crease



Description	Dimensions (W × H)	Colour	Unit	Order no.
<b>70 Shore A</b>				
Cushion Crease	16.5 × 6.4 mm	Yellow	1 Coil / 50 m	BPB00142
Cushion Crease	17.0 × 7.0 mm	Green	1 Coil / 50 m	BPB00143
Cushion Crease	17.5 × 7.5 mm	Gold	1 Coil / 50 m	BPB00153
Cushion Crease	18.0 × 8.0 mm	Blue	1 Coil / 50 m	BPB00144
Cushion Crease	18.5 × 8.5 mm	Purple	1 Coil / 50 m	BPB00151
Cushion Crease	20.0 × 9.0 mm	Orange	1 Coil / 50 m	BPB00145
Cushion Crease	20.0 × 9.7 mm	White	1 Coil / 50 m	BPB00152
Cushion Crease	20.0 × 10.0 mm	Pink	1 Coil / 50 m	BPB00154


Note: The indicated Shore A hardnesses can vary slightly.

### Cushion Crease EasyFix




Description	Dimensions (W × H)	Colour	Unit	Order no.
<b>70 Shore A</b>				
Cushion Crease EasyFix	16.5 × 6.4 mm	Yellow	1 Roll / 20 m	BPB01242
Cushion Crease EasyFix	17.0 × 7.0 mm	Green	1 Roll / 20 m	BPB01243
Cushion Crease EasyFix	17.5 × 7.5 mm	Gold	1 Roll / 20 m	BPB01253
Cushion Crease EasyFix	18.0 × 8.0 mm	Blue	1 Roll / 20 m	BPB01244
Cushion Crease EasyFix	18.5 × 8.5 mm	Purple	1 Roll / 18 m	BPB01251
Cushion Crease EasyFix	20.0 × 9.0 mm	Orange	1 Roll / 16 m	BPB01245
Cushion Crease EasyFix	20.0 × 9.7 mm	White	1 Roll / 15 m	BPB01252
Cushion Crease EasyFix	20.0 × 10.0 mm	Pink	1 Roll / 15 m	BPB01254

Note: The indicated Shore A hardnesses can vary slightly.




## CITO Cushion Crease Plus – creasing profile



CITO Cushion Crease Plus is mainly used for double and triple flute paper. It is placed on the creasing rule to support compression in a larger area of the corrugated board and to optimise folding behaviour, particularly along the corrugations.


### Cushion Crease Plus



Description	Dimensions (W × H)	Colour	Unit	Order no.
<b>70 Shore A</b>				
Cushion Crease Plus	21 × 7.0 mm	Green	1 Coil / 50 m	BPB00171
Cushion Crease Plus	21 × 7.5 mm	Gold	1 Coil / 50 m	BPB00172
Cushion Crease Plus	21 × 8.0 mm	Blue	1 Coil / 50 m	BPB00173
Cushion Crease Plus	21 × 8.5 mm	Purple	1 Coil / 50 m	BPB00174
Cushion Crease Plus	21 × 9.0 mm	Orange	1 Coil / 50 m	BPB00175
Cushion Crease Plus	21 × 9.5 mm	White	1 Coil / 50 m	BPB00176
Cushion Crease Plus	21 × 11.0 mm	Beige	1 Coil / 50 m	BPB00180
Cushion Crease Plus	21 × 12.0 mm	Brown	1 Coil / 35 m	BPB00182

Note: The indicated Shore A hardnesses can vary slightly.

### Cushion Crease Plus EasyFix



Description	Dimensions (W × H)	Colour	Unit	Order no.
<b>70 Shore A</b>				
Cushion Crease Plus EasyFix	21 × 9.0 mm	Orange	1 Roll / 16 m	BPB01275
Cushion Crease Plus EasyFix	21 × 9.5 mm	White	1 Roll / 14 m	BPB01276
Cushion Crease Plus EasyFix	21 × 11.0 mm	Beige	1 Roll / 12 m	BPB01280

Note: The indicated Shore A hardnesses can vary slightly.



CITO Cushion Crease and CITO Cushion Crease Plus are suitable for use on **flatbed and rotary dies**.

The creasing profile is characterised by its **excellent height and wear resistance**.

Due to the **colour-coding**, the quality control process can be carried out by simply visually inspecting the cutting die.

CITO Cushion Crease and CITO Cushion Crease Plus are available as a **self-adhesive EasyFix** version on **user-friendly CMR reels**. Thanks to the EasyFix technology, the ejection material can be precisely positioned on the cutting die and can be removed without leaving any residues. The surface of your cutting die remains undamaged. Additional stapling is not necessary.



#### CMR Holder for CITO Creasing Matrix – for the efficient storage and conversion of creasing systems on reels

The compact CMR Holder for CMR Discs is characterised by its stable design and can hold up to four CMR reels.

Advantages at a glance:

- Efficient storage and conversion of CMR reels
- Compact and space-saving design
- Stable construction
- Suitable for up to four CMR reels
- CMR reels are individually removable

Optimally suitable for:

- CITO ULTIMATE CMR
- CITO DYNAMIC CMR
- CITO BASICplus CMR
- CITO Cushion Crease EasyFix and CITO Cushion Crease Plus EasyFix

Description	Unit	Order no.
CITO CMR Holder, 4-Slot	1 piece	BDA00002

Dimensions: Depth: 460 mm, Width: 380 mm, Height: 410 mm

Installation dimensions incl. CMR disc: Depth: 530 mm, Width: 380 mm, Height: 655 mm



Further information  
on CITO Cushion Crease:  
[www.cito.de/cushion-crease](http://www.cito.de/cushion-crease)